




FILTER TECHNOLOGY

MOLECULAR BIOLOGY PRODUCT COLLECTION





FILTER TECHNOLOGY



Molecular Biology

Index


Blotting Overview	2
--------------------------------	----------

Membranes

Pure Nitrocellulose	4
Supported Nitrocellulose	5
Polyvinylidene Fluoride PVDF	6


ECL Substrates

Lightwave	7
Lightwave Plus	10
Lightwave Max	12




MOLECULAR BIOLOGY ANALYSIS

Molecular analysis studies subcellular components such as proteins and nucleic acids (DNA, RNA). These molecules can be detected by various blotting techniques. The sample of interest is separated according to size by electrophoresis through a gel. Molecules from the sample are transferred and bound to a microporous membrane. Then, specific molecules of interest are detected using another molecule which specifically binds to the molecule of interest and can be detected by color, light or radioactivity.



Western Blot

Western blotting is a common and important technique used in molecular biology. It is used to detect a specific protein or protein fragment from a complex mixture such as a cell lysate, tissue extract, blood or serum sample or culture supernatants. The complex mixture is separated according to size by gel electrophoresis and then transferred to a membrane. A protein of specific interest is immunodetected using primary and secondary antibodies.



Western Blot Application Examples:

- Protein expression and modification studies, may be quantitative;
- Amino acid analysis;
- Diagnostics development;
- Medical diagnosis such as for HIV and Lyme disease.

WESTERN BLOTTING PROTOCOL

Electrophoretic separation of proteins

Separation into polyacrylamide gel according to molecular weight. In order to separate the proteins of lower molecular weight, use of more concentrated gel is required.


Transfer of proteins

Transfer from gel onto membrane followed by:

- Blocking;
- Applying a primary antibody specific for your protein of interest;
- Applying secondary antibody that will recognize the primary antibody.

Role of protein binding

Set up for transfer



TYPICALLY TWO TYPES OF MEMBRANES:


- Nitrocellulose (hydrophilic)
- Polyvinylidene fluoride (hydrophobic)

Protein to membrane binding interactions: hydrophobic, electrostatic, dipolar



Detection of proteins

Proteins can be detected by immunodetection methods which use enzyme conjugated/labeled secondary antibodies. When the enzyme substrate is added, a product is formed. This product can be detected by fluorescence, colorimetrically, or by chemiluminescence. Enhanced chemiluminescence (ECL) produces light as a by-product when the substrate is catalyzed by the enzyme. This light is then captured on X-ray film or by a digital imaging system.



No separation of proteins

The protein is spotted on to a nitrocellulose or PVDF membrane, for example using a low-volume pipette.



Serial dilution of purified proteins




Pure Nitrocellulose



Pure Nitrocellulose is the membrane of choice for all protein or immunoblotting application. The most common used membrane for western blotting techniques. Supplied in various porosity and format


Features and Benefit

- ◆ High resolution
- ◆ Low Background, easily blocked
- ◆ Wets out naturally
- ◆ Compatible with all detection system



Ordering information

	Dimensions (mm) Packaging	70x84 mm 10/pk	100x100 mm 10/pk	150x150 mm 5/pk	200x200 mm 25/pk	200x3000 mm 1/pk	300x3000 mm 1/pk
Pore sizes	0.22 µm	1213991	1213999	1215463	1215392	1215469	1215458
	0.45 µm	1213888	1213314	1215476	1221976	1215483	1215471



Supported Nitrocellulose




Supported Nitrocellulose Transfer Membrane combines the binding characteristics of nitrocellulose membrane with the strength of nylon membrane.

Supplied in various porosity and format

Features and Benefit

- ▲ Supported for procedures requiring rigorous handling
- ▲ Strong will not curl, bend or crack after baking
- ▲ High sensitivities, low backgrounds
- ▲ Multiple reprobings
- ▲ BSA binding capacity up to



All lanes : Anti-Furin antibody [EPR14674] (ab183495) at 1/5000 dilution

Lane 1 : HepG2 whole cell lysate

Lane 2 : HeLa whole cell lysate

Lane 3 : U87-MG whole cell lysate

Lane 4 : Caco-2 whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size : 87 kDa

Ordering information

Dimensions (mm) Packaging	70x84 mm 10/pk	100x100 mm 10/pk	150x150 mm 5/pk	200x200 mm 5 /pk	200x3000 mm 1/pk	300x3000 mm 1/pk
0.22 µm		1214560	1212669	1212689	1212690	1212632
0.45 µm	1214978	1213943	1212596	1212597	1212602	1212590



Polyvinylidene Fluoride PVDF



PVDF is a naturally hydrophobic unsupported transfer membrane.

It has a high binding capacity and low backgrounds

Supplied in various porosity and format


Features and Benefit

- ▲ Superior Strength: Can withstand aggressive handling or be used with automated equipment without breaking or tearing
- ▲ Low extractable: Ensures tests will be clean with consistent results
- ▲ Exceptional sensitivity: Detects low-level components
- ▲ Hydrophobic: For high protein binding
- ▲ Lot to lot consistency: Quality checks ensure consistent binding for dependable results every time
- ▲ High binding capacity: Binds a wide range of fragment sizes
- ▲ High range of chemical Resistant to most commonly used chemicals compatible with chemically aggressive solvents

Ordering information

Dimensions (mm) Packaging	70x84 mm 10/pk	100x100 mm 10/pk	150x150 mm 5/pk	200x200 mm 5/pk	200x3000 mm 1/pk	300x3000 mm 1/pk
0.22 µm	1214588		1215037	1215032	1214726	1214429
0.45 µm	1213992	1212644	1212636	1212637	1212783	1212639

Pore sizes



LIGHTwave™

ECL SUBSTRATES FOR WESTERN BLOTTING

About us

The GVS Group is one of the world's leading manufacturers of filters and components for applications in the Healthcare, Life Sciences, Automotive, Appliance, Safety, and Commercial & Industrial Filtration.

The Group's clear strategy towards internationalization, has led to the opening of 12 production facilities located in Italy, UK, Brazil, the United States, China and Romania, as well as offices in Russia, Turkey, Argentina, Japan, Korea. GVS currently have a workforce of over 2,700 people globally.


For 40 years, GVS has focused on innovation in its products range and production processes, constantly improving its development capacity to provide the best service and support for its clients

We offer a full range of branded products through a global network of dealers and distributors. We also make available all these capabilities on an OEM basis by working closely with companies around the world to provide state of the art materials solutions and/or turn-key final product solutions used in critical applications for the pharmaceutical, medical device, diagnostic, food & beverage and environmental monitoring markets.

All GVS substrates are protected by **US7803573**, **EP1962095**, **US7855287**, **EP1950207**, **US2012009603 (A1)**, **CA2742025**, **EP2405016**, foreign equivalents and pending patents.

LightWave™ detection reagents are non-isotopic, luminol-based chemiluminescence substrate, designed for the chemiluminescent detection of immobilized proteins and immobilized nucleic acids conjugated with horseradish peroxidase (HRP).

LightWave™ is intended for research use only, and shall not be used in any clinical procedures, or for diagnostic purposes.



Storage/expiry

One year at room temperature (max. 25°C).


LightWave™ product line

Product	LightWave™	LightWave™ Plus	LightWave™ Max
Signal intensity	Medium	High	Ultra High
Signal duration	Medium	Extended	Short
Protein abundance	High	Medium	Ultra-low

GVS LIGHTWAVE SUBSTRATES

Overview

HPR - Antibody dilutions



LightWave - Low picogram detection level
 LightWave Plus - Mid femtogram detection level
 LightWave Max - low femtogram detection level

Product	Suggested antibody dilutions	
LightWave™	Primary Ab	1:500 - 1:5,000
	Secondary Ab	1:20,000 - 1:100,000
LightWave™ Plus	Primary Ab	1:1000 - 1:15,000
	Secondary Ab	1:25,000 - 1:150,000
LightWave™ Max	Primary Ab	1:5000 - 1:100,000
	Secondary Ab	1:100,000 - 1:500,000

Quick start protocol

- ◆ Perform electrophoresis, membrane transfer and antibody incubation and washes
- ◆ Prepare Lightwave™ ECL substrate by mixing equal volumes of the two solutions
- ◆ Apply Lightwave™ chemiluminescent substrate to the membrane (1 mL per 10 cm³ of the membrane), incubate 2 minutes with the substrate
- ◆ Expose the substrate-treated membrane using a chemiluminescence imager or X-ray film

Product	Competitors
LIGHTwave™	PIERCE™ ECL PLUS - THERMO SCIENTIFIC™ IMMOBILION® CLASSICO - MILLIPORE™ WESTERN LIGHTNING™ PLUS - PERKINELMER WESTERNBRIGHT™ ECL - ADVANSTA
LIGHTwave™ Plus	CLARITY™ - BIORAD SUPERSIGNAL™ WEST DURA - THERMO SCIENTIFIC™ AMERSHAM™ ECL PRIME™ - GE HEALTHCARE SUPERSIGNAL™ WEST PICO PLUS - THERMO SCIENTIFIC™ IMMOBILION® CRESCENDO - MILLIPORE™ WESTERNBRIGHT™ QUANTUM™ - ADVANSTA
LIGHTwave™ Max	CLARITY MAX™ - BIORAD SUPERSIGNAL™ WEST FEMTO - THERMO SCIENTIFIC™ AMERSHAM™ ECL SELECT™ - GE HEALTHCARE WESTERNBRIGHT™ SIRIUS™ - ADVANSTA WESTERN LIGHTNING™ ULTRA - PERKINELMER

GVS Lightwave



LIGHTwave™

Competitor Pico

Competitor Classico


Competitor ECL

Features

- ◆ Low picogram detection
- ◆ Ideal for routine analysis
- ◆ Working solution stable for at least three days
- ◆ The best entry level ECL substrate on the market
- ◆ Signal duration 5 hours
- ◆ Stable for 1 year at RT

Code	Description
LW0001	LightWave™ Western blotting substrate 10 mL
LW0002	LightWave™ Western blotting substrate 250 mL

GVS Lightwave Plus




Plus
LIGHTwave™

Competitor B

Features

- ◆ Mid femtogram detection
- ◆ Extended signal duration
- ◆ High range flexibility
- ◆ Working solution stable for at least three days
- ◆ The perfect ECL formulation combining great sensitivity and long signal duration
- ◆ Signal duration 25 hours
- ◆ Stable for 1 year at RT




Plus
LIGHTwave™

Competitor

Western blotting detection of HDAC-1 on HeLa cells

HeLa cell lysate from 5 to 0,078 µg
 Ab 1° Rabbit anti HDAC1 1:5000
 Ab 2° Goat anti rabbit 1:50000
 Exposure time: 3 minutes
 Imager: LAS4000 (GEHC)


Signal to noise ratio



GVS LightWave Plus vs Competitor Signal duration

Signal duration


LightWave™ Plus provides an extremely extended signal duration when compared to most mid-level range ECL substrates. The HDAC-1 signal intensity variation over time was analyzed using **LightWave™ Plus** and its competitors (Figure 3).



Code	Description
LW0003	LightWave™ Plus Western Blotting Substrate 10 ml
LW0004	LightWave™ Plus Western Blotting Substrate 250 ml



GVS Lightwave Max



Max LIGHTwave™

Competitor Femto

Competitor ECL Select

Features


- ◆ Low femtogram detection
- ◆ Low antibody consumption to save your money
- ◆ Working solution stable for at least three days
- ◆ The ECL substrate with the highest signal on the market
- ◆ Signal duration 8 hours
- ◆ Stable for 1 year at RT



Figure 2. Low background for high sensitive detection with LightWave™ Max.

A) Western blotting detection of HDAC-1 on HeLa cell lysate with LightWave™ Max compared to SuperSignal™ West Femto-Thermo Scientific™. Triplicate blots for each substrate containing 2-fold dilutions of HeLa whole cell lysate were incubated with primary antibody (Rabbit-anti Human HDAC-1) 1:15000 and secondary antibody (Goat anti Rabbit-HRP) 1: 300000 and were simultaneously imaged for 120 seconds with ImageQuant™ LAS 4000 (GE Healthcare).

B) Signal-to-noise ratio (S/N) analysis. LightWave™ Max displays the best combination of sensitivity and signal with low background.



Detection level: Low-femtogram




Figure 1. Western blotting detection of HDAC-1 on HeLa cell lysate with LightWave™ Max and other chemiluminescent reagents in the same sensitivity range.

Signal duration




Figure 3. Signal duration of LightWave™ Max compared to SuperSignal™ West Femto-Thermo Scientific™.


Quadruplicate blots for each substrate containing 2-fold dilutions of HeLa whole cell lysate were incubated with primary antibody (Rabbit-anti Human HDAC-1) 1:15000 and secondary antibody (Goat anti Rabbit-HRP) 1: 300000 and were simultaneously imaged with ImageQuant™ LAS 4000 (GE Healthcare) at time points up to 11 hours post substrate addition.

Code	Description
LW0005	LightWave™ Max Western Blotting High Sensitive Substrate 10 ml
LW0006	LightWave™ Max Western Blotting High Sensitive Substrate 100 ml



WORLDWIDE
DISTRIBUTION CENTERS

FILTER TECHNOLOGY



EUROPE

Italy Office
Headquarters
GVS S.p.A.
Via Roma 50
40069 Zola Predosa (BO) - Italy
tel. +39 051 6176311
fax +39 051 6176200
gvs@gvs.com

United Kingdom
GVS Filter Technology UK Ltd.
NFC House
Vickers Industrial Estate
Mellishaw Lane, Morecambe
Lancashire LA3 3EN
tel. +44 (0) 1524 847600
lifesciences.uk@gvs.com

Russia
GVS Russia LLC
Profsoyuznaya Street, 25-A, office 102
117418, Moscow
Russian Federation (Russia)
tel. +7 495 0045077
lifesciences.ru@gvs.com

Romania
GVS Microfiltrazione srl
Str. Principala n. 320 et. 1 – Ciorani de Jos
JUD . PRAHOVA - CIORANI ROMÂNIA
Tel. +40 244 463044
lifesciences.ro@gvs.com

Turkey
GVS Türkiye
Cevizli mah. Zuhul cad. Ritim Istanbul
no: 44 A-1 Blok D.371 Maltepe / Istanbul
tel. +90 216 504 47 67
lifesciences.tr@gvs.com

ASIA

China
GVS Technology (Suzhou) Co., Ltd.
Fengqiao Civil-Run Sci-Tech Park,
602 Changjiang Road,S.N.D.
Suzhou, China 215129
tel. +86 512 6661 9880
fax: +86 512 6661 9882
lifesciences.cn@gvs.com

Japan
GVS Japan K.K.
KKD Building 4F, 7-10-12 Nishishinjuku
Shinjuku-ku, Tokyo 160-0023 Japan
tel. +81 3 5937 1447
fax +81 3 5937 1448
lifesciences.jp@gvs.com

Korea
GVS Korea Ltd #315 Bricks Tower
368 Gyungchun-ro(Gaun-dong),
Namyangju-si, Gyunggi-do,
Tel: +82 31 563 9873
Fax: +82 31 563 9874
lifesciences.kr@gvs.com

AMERICA

U.S.A.
GVS North America, Inc.
63 Community Drive
Sanford, ME 04073 - USA
tel. +1 866 7361250
lifesciences.us@gvs.com

Mexico
GVS de México
Universal No. 550, Vynmsa Aeropuerto Apodaca
Industrial Park, Ciudad Apodaca, Nuevo León,
C.P. 66626 México
tel. +52 81 2282 9003
lifesciences.mx@gvs.com

Brazil
GVS do Brasil Ltda.
Rodovia Conego Cyriaco Scaranello Pires 251
Jd. Progresso, CEP 13190-000
Monte Mor (SP) - Brasil
tel. +55 19 38797200
fax +55 19 38797251
lifesciences.br@gvs.com

Argentina
GVS Argentina S.A.
Francisco Acuña de Figueroa
719 Piso:11 Of: 57
1416 Buenos Aires - Argentina
tel. +54 11 49889041
fax +54 11 49889042
lifesciences.ar@gvs.com